

Clock Configuration (Forensic Models Only)

To configure the internal clock module, while in Flash Drive Mode, open the file named TIME_OLD.TXT. Update the file using the following format and then *Save As*: TIME.TXT

```
Year=2009  
Month=04  
Day=01  
Hour=12  
Minute=34  
Second=56  
Format=PM
```

Values must match the options below:

- Year sets the clock year value. (2000 to 2099)
- Month sets the clock month value. (01 to 12)
- Day sets the clock day value. (01 to 31)
- Hour sets the clock hour value. (01 to 12) or (00 to 23)
- Minute sets the clock minute value. (00 to 59)
- Second sets the clock second value. (00 to 59)
- Format sets the time format. (AM, PM, or 24) If AM is chosen, the 12-hour format is selected and the specified hour is treated as before noon. If PM is chosen, the 12-hour format is selected and the specified hour is treated as after noon. If 24 is chosen, the 24-hour format is selected and the specified hour is treated as 24-hour format.

Note: The time and date will be updated when the keylogger is reconnected. For precise time setting, it is a good idea to set the time a few minutes in advance, and restart the keylogger at the specified time.